

# CLINICAL EVIDENCE FOR ROBOTICS IN EP

## DATA TRANSPARENCY

Because safety and better outcomes for patients are your first priorities



ROBOTIC VS. MANUAL







>50 PATIENTS



2012-2020

We believe data should be presented clearly and comprehensively. That's why we present data from all known publications describing **robotic** vs. manual cardiac catheter ablations in head-to-head comparisons with greater than 50 patients since 2012.



Stereotaxis robotic technology has been described in more than 400 scientific publications. Review the above data in detail and discover the full, searchable database of publications at [www.RoboticEP.com/clinical-data](http://www.RoboticEP.com/clinical-data)



# KEY CLINICAL EVIDENCE FOR ROBOTICS IN EP

## ATRIAL FIBRILLATION

- Ghadban R et al. .2020. "Radiofrequency atrial fibrillation with irrigated tip catheter using remote magnetic navigation compared with conventional manual method." *Journal of Interventional Cardiac Electrophysiology* 95-102.
  - MAN & RMN are comparable in terms of AF Ablation success rates
  - RMN has lower complication rates and fluoroscopy times, but longer procedural times.
- Noten AME et al. 2019. "Robotic Navigation Shows Superior Improvement in Efficiency for Atrial Fibrillation Ablation." *Journal of Atrial Fibrillation*
  - RMN showed greater generational improvement in procedure times compared to CRYO and MCN
- Luo Q et al. "Utilization of steerable sheath improves the efficiency of atrial fibrillation ablation guided by robotic magnetic navigation compared with fixed-curve sheath." *Clin Cardiol.* 2022 Feb 23. doi: 10.1002/clc.23801. Epub ahead of print. PMID: 35195273.
  - Steerable sheaths reduce procedure times versus fixed sheaths:  $90.4 \pm 20.7$  versus  $111.9 \pm 25.2$  min,  $p < .001$
  - Utilization of steerable sheath technology can improve the efficiency of AF ablation guided by RMN, primarily by reducing the total procedure and RF delivery times without compromising safety. Procedure time  $111.9 \pm 25.2$  vs.  $90.4 \pm 20.7$  min,  $p < .001$

## VENTRICULAR ARRHYTHMIAS

- Blandino A, et al. "Outcomes of manual versus remote magnetic navigation for catheter ablation of ventricular tachycardia: a systematic review and updated meta-analysis". *Pacing Clin Electrophysiol.* 2021 Jun;44(6):1102-1114. doi: 10.1111/pace.14231. Epub 2021 Apr 24. PMID: 33825206.
  - RMN was superior to MAN guided CF ablations in terms of acute success, complications, and fluoro times
- Hendriks et al. "Safety & Clinical Outcome of Catheter Ablation of Ventricular Arrhythmias Using Contact Force Sensing, Consecutive Case Series" *J Cardiovasc Electrophysiol.* 2015 Jul
  - RMN was superior with regard to acute success, reduction of major complications, and recurrence rate.

## SVT & CONGENITAL

- Roy K et al. 2016 "RMN for catheter ablation in patients with CHD: A Review" *Journal of Cardiovascular Electrophysiology.* Vol 6
  - Editorial on advantages of RMN in the setting of ACHD
- Noten AME et al. 2021. "RMN shows superior long-term outcomes in pediatric atrioventricular (nodal) tachycardia ablation compared to manual radiofrequency and cryoablation" *International Journal of Cardiology Heart and Vascular.* Vol 3
  - RMN had significantly lower recurrence rates compared to MAN and CRYO at mean follow-up of  $5.5 \pm 2.9$  years (AVRT: 4.3% versus 15.6% versus 54.5%,  $P < 0.001$ ; AVNRT: 7.7% versus 8.3% versus 35.7%,  $p = 0.008$ ; for RMN versus MAN versus CRYO respectively).

## CLINICAL SCIENCE

- Bhaskaran A et al. "Magnetic guidance versus manual control: comparison of radiofrequency lesion dimensions and evaluation of the effect of heart wall motion in a myocardial phantom." *J Interv Card Electrophysiol.* 2015;44(1):1-8.
  - In vitro simulations show that robotic catheters are more stable on simulated wall motion, leading to more focal lesions
- Bun S-S et al. 2017. "Radiofrequency catheter ablation of atrial fibrillation: Electrical modification suggesting transmural is faster achieved with remote magnetic catheter in comparison with contact force use." *Journal of Cardiovascular Electrophysiology* 1-9.
  - Faster EGM modification is achieved with RMN due to increased catheter stability
- Luo Q et al. "Utilization of steerable sheath improves the efficiency of atrial fibrillation ablation guided by robotic magnetic navigation compared with fixed-curve sheath." *Clin Cardiol.* 2022 Feb 23. doi: 10.1002/clc.23801. Epub ahead of print. PMID: 35195273.
  - Signal energy attenuation and impedance drop were evaluated as lesion dimension surrogates and signal fragmentation and shrinkage were estimated as lesion quality surrogates.
  - RMN creates more homogeneous lesion sets. MCN lesions are larger but less homogeneous relative to RMN.

